


**Light Collection Efficiency
in Thin Strip Plastic
Scintillator for the Study of
ISGMR in Unstable Nuclei**




Giant Resonances

- Collective excitations of the nuclei
- Discovered in the 1940s while bombarding nuclei with gamma rays
- Monopole resonance is a spherical oscillation
- Isoscalar – neutrons and protons move in phase with one another
- Isovector – neutrons and protons move out of phase with one another




Isoscalar Giant Monopole Resonance (ISGMR)

- ISGMR: the breathing mode
 - From studying the Isoscalar Giant Monopole Resonance (ISGMR), E_{GMR} can be determined
 - E_{GMR} can be used to find K_A ,
 - K_A can be used to find K_{NM}
- K_{nm} is a quantity which describes a ground state property of nuclear matter
- Uses include:
 - Determining the nuclear equation of state
 - Neutron Stars
 - Supernova collapse
 - Heavy ion collisions



$$E_{GMR} = \sqrt{\frac{\hbar^2 A K_A}{m \langle r^2 \rangle}}$$




Unstable Nuclei



ISGMR in unstable nuclei:

- Determine the dependence of the nuclear incompressibility on $(N-Z)/A$
- Cannot successfully make a target out of unstable nuclei
- Study the inverse reaction

Normal Kinematics:



Inverse Kinematics:



Experimental Setup


Decay detector:

- Detects decay by proton and alpha emission
- Particle type
- Angle
- Energy

Gas detector:

- Detects heavy ions
- Energy
- Angle

MDM




Measuring Protons and Alphas from Decay of ISGMR

- Detector:
 - Layer of horizontal 1 mm thick scintillator strips
 - Layer of vertical 1 mm thick scintillator strips
 - 5 block scintillators
- Scintillator strips:
 - Angle
 - Energy loss
- Block scintillators:
 - Total energy
- Angular resolution $\sim 4^\circ$



Scintillators


- Scintillator:
 - absorbs energy
 - emits light
- Particle:
 - Excites electrons
 - Photons are emitted as electrons return to low energy states.
- Photons are translated into an electrical signal by a photomultiplier which is coupled to the scintillator.




Calculation of light output in $\Delta E1$ and $\Delta E2$ for protons and alpha ions

Photomultiplier

- Photomultiplier:
 - Absorbs light
 - Releases Electrons via photoelectric effect at the photocathode
 - The cathode, dynodes, and anodes create a potential “ladder”
 - Electrons travel from the photocathode to the first dynode and excite more electrons in the dynode
 - Excited electrons leave the dynode and travel to the next dynode in sequence and repeats the process
 - At the anode all the electrons are collected and then amplify to create a readable current
- Increasing light collection efficiency:
 - Wrap the scintillators in reflective material




Scintillator Wrapping


- Wrapping:
 - Light tight
 - Reflective
 - Reflective material:
 - Ensures the maximum number of photons enter the photomultiplier.
 - Optimize the improved signal vs. the energy lost by the particle.
- 

Types of Wrapping


- Aluminum foil
- Aluminized Mylar
- 3M Tyvek
- reflective paint



Scintillator Testing Setup




^{241}Am Spectrum




Typical ^{241}Am Spectrum Obtained with an ORTEC Partially-Depleted Detector.


^{241}Am Spectrum



Spectrum Shift with Distance



Data




For 0 layer $A_0 = 675 \pm 7$
 $B_0 = 37.6 \pm 0.8$

For 1 layer $A_1 = 661 \pm 13$
 $B_1 = 30.0 \pm 1.5$


Testing

- Light proof box
- Scintillator to PMT coupling
- Calibrating electronics



Continuing Work

- Test various wrapping to determine the most effective type to use in the detector
- Assemble, test and use detector to study ISGMR in various unstable nuclei



Acknowledgements

Thank you to Dr. Youngblood, Dr. Lui and Dr. Krishichayan for their guidance during this project, the Cyclotron Institute and Dr. Yennello for hosting the REU program, NSF for funding this opportunity, and a special thanks to Jon Button for all of the extra time he put into helping me this summer.

